



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human JUN Polyclonal Antibody, Rabbit

Product Code: KB5556GNP

Gene: JUN (Transcription factor AP-1, Activator protein 1, AP1, Proto-oncogene c-jun, V-jun avian sarcoma virus 17 oncogene homolog, p39)

Immunogen: GX5102 (GST-fusion protein, 163 amino acids)
KMETTFYDDALNASFLPSESGPYGYSNPKILKQSMTLNLADPVGSLKP
HLRAKNSDLLTSPDVGLLKLASPELERLIQSSNGHITTTPTPTQFLCPKN
VTDEQEGFAEGFVRALAEHSQNTLPSVTSAAPVNGAGMVAPAVASV
AGGSGSGGFSASLHSE

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human JUN protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619