



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human HMG20A Polyclonal Antibody, Rabbit

Product Code: KB5656GNP

Gene: HMG20A (High mobility group protein 20A, HMG box-containing protein 20A, HMG domain-containing protein HMGX1)

Immunogen: GX5277 (GST-fusion protein, 172 amino acids)
YKVFSRKTQDRQKGKSHRQDAARQATHDHEKETEVKERSVFDIPIFTE
EFLNHSKAREAELRQLRKSMEFEERNAALQKHVESMRTAVEKLEVDV
IQERSRNTVLQQHLETLRQVLTSSFASMPPLPGSGETPTVDTIDSYMNRL
HSIILANPQDNENFIATVREVVNRLDR

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human HMG20A protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619