



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human TRIM3 Polyclonal Antibody, Rabbit

Product Code: KB5705GNP

Gene: TRIM3 (Tripartite motif-containing protein 3, RING finger protein 22, Brain-expressed RING finger protein, RING finger protein 97)

Immunogen: GX5274 (GST-fusion protein, 176 amino acids)
KIGAGRLMGPKGVAVDRNGHIIVVDNKSCCVFTFQPNGKLVGRFGGRG
ATDRHFAGPHFVAVNNKNEIVTDFHNHVKVYSADGEFLFKFGSHGE
GNGQFNAPTGVAVDSNGNIIVADWGNRIQVFDSSGSFLSYINTSAEPL
YGPQGLALTSDGHVVADAGNHCFKAYRYLQ

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human TRIM3 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619