



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human PSMC5 Polyclonal Antibody, Rabbit

Product Code: KB6588GNP

Gene: PSMC5 (26S protease regulatory subunit 8, Proteasome 26S subunit ATPase 5, Proteasome subunit p45, p45/SUG, Thyroid hormone receptor-interacting protein 1, TRIP1)

Immunogen: GX5159 (GST-fusion protein, 163 amino acids)
AGSGLRQYYLSKIEELQLIVNDKSQNLRLQAQRNELNAKVRLREELQ
LLQEQGSYVGEVVRAMDKKKVLVKVHPEGKFVVDVDKNIDINDVTPNC
RVALRNDSTLHKILPNKVDPLVSLMMVEKVPDSTYEMIGGLDKQIKEIK
EVIELPVKHPELFEAL

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human PSMC5 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619