



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human TNRC4 Polyclonal Antibody, Rabbit

Product Code: KB7008GNPAF

Gene: TNRC4 (CUG-BP- and ETR-3-like factor 3, CELF-3, Bruno-like protein 1, RNA-binding protein BRUNOL-1, ELAV-type RNA-binding protein 1, ETR-1, Trinucleotide repeat-containing gene 4 protein, Expanded repeat domain protein CAG/CTG 4, CAG repeat protein 4)

Immunogen: GX5067 (GST-fusion protein, 167 amino acids)
MKEPDAIKLFVGQIPRHLEEKDLKPIFEQFGRIFELTVIKDKYTGLHKGCA
FLTYCARDSALKAQSALHEQKTLPGMNRPIQVKPADSESERGEDRKLFV
GMLGKQQTDEDVRKMFEPFGTIDECTVLRGPDGTSKGCAFVKFQTHA
EAQAAINTLHSSRTLPGASSS

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human TNRC4 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619