



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ANP32A Polyclonal Antibody, Rabbit

Product Code: KB7264GNP

Gene: ANP32A (Acidic leucine-rich nuclear phosphoprotein 32 family member A, Potent heat-stable protein phosphatase 2A inhibitor I1PP2A, Acidic nuclear phosphoprotein pp32, Leucine-rich acidic nuclear protein, Lanp, Putative HLA-DR-associated protein I, PHAPI, Mapmodulin)

Immunogen: GX5388 (GST-fusion protein, 153 amino acids)
LELRNRTPSDVKELVLDNSRSNEGKLEGLTDEFEELEFLSTINVGLTSIA
NLPKLNKLLKLELSDNRVSGGLEVLAEKCPNLTHLNLSGNKIKDLSTIEP
LKKLENLKSLEDFNCEVTNLNDYRENVFKLLPQLTYLDGYDRDDKEAPD
SDAE

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human ANP32A protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619