



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human SAP18 Polyclonal Antibody, Rabbit

Product Code: KB8109GNP

Gene: SAP18 (Histone deacetylase complex subunit SAP18, Sin3-associated polypeptide, 18 kDa, Sin3-associated polypeptide p18, Cell growth-inhibiting protein 38, 2HOR0202)

Immunogen: GX5395 (GST-fusion protein, 151 amino acids)
MAVESRVTQEEIKKEPEKPIDREKTCPLLLRVFTTNNGRHHRMDEF SR
GNVPSSELQIYTWMDATLKELTSLVKEVYPEARKKGT HFNFAIVFTDVK
RPGYRVKEIGSTMSGRKGTDDSM TLQSQKFQIGDYLDIAITPPNRAPPP
SGRMR

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human SAP18 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619