



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ZSCAN18 Polyclonal Antibody, Rabbit

Product Code: KB8132GNP

Gene: ZSCAN18 (Zinc finger and SCAN domain-containing protein 18, Zinc finger protein 447)

Immunogen: GX5291 (GST-fusion protein, 96 amino acids)
GILPDKVRPWVVAQYPESCKKAASLVEGLADVLEEPGMMLLGSPAGSSSI
LSDGVYERHMDPLLLPGELASPSQALGAGEIPAPSETRECPASLAAS

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human ZSCAN18 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619