



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human OPTN Polyclonal Antibody, Rabbit

Product Code: KB9422GNP

Gene: OPTN (Optineurin, Optic neuropathy-inducing protein, E3-14.7K-interacting protein, FIP-2, Huntingtin-interacting protein L, Huntingtin yeast partner L, NEMO-related protein, Transcription factor IIIA-interacting protein, TFIIIA-IntP)

Immunogen: GX5384 (GST-fusion protein, 176 amino acids)
EHNNALKTIEELTRKESEKVDRAVLKELSEKLELAEKALASKQLQMDDEM
KQTIKQEEDLETMTILRAQMEVYCSDFAERAAREKIHEEKEQLALQL
AVLLKENDAFEDGGRQSLMEMQSRHGARTSDSDQQAYLVQRGAEDR
DWRQQRNIPHIHSCPKCGEVLDPDIDTLQIHVMD

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human OPTN protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619