



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human GTF2F1 Polyclonal Antibody, Rabbit

Product Code: KB9482GNP

Gene: GTF2F1 (General transcription factor IIF subunit 1, Transcription initiation factor IIF subunit alpha, TFIIIF-alpha, Transcription initiation factor RAP74, General transcription factor IIF polypeptide 1, 74 kDa subunit protein)

Immunogen: GX5554 (GST-fusion protein, 174 amino acids)
ESDSSEESDIDSEASSALFMAKKKTPPKRERKPSGGSSRGNSRPGTPS
AEGGSTSSTLRAAASKLEQGKRVSEMPAAKRLRLDTGPQSLSGKSTP
QPPSGKTTTPNSGDVQVTEDAVRRYLTRKPMTTKDLLKKFQTKKTGLSS
EQTNVLAQILKRLNPERKMINDKMHFSLKE

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human GTF2F1 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619