



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human KLHL25 Polyclonal Antibody, Rabbit

Product Code: KB9611GNP

Gene: KLHL25 (Ectoderm-neural cortex protein 2, ENC-2, Kelch-like protein 25)

Immunogen: GX5438 (GST-fusion protein, 164 amino acids)
MSVSVHETRKSRSTGSMNVTLFHKASHPDCVLAHLNLTNRKHCMTD
VTLWAGDRAFPCHRAVLAASSRYFEAMFSHGLRESRDDTVNFQDNLH
PEVLELLLDFAVSSRIANEENAESLLEAGDMLQFHDVRDAAAEFLEKNL
FPSNCLGMMLLSDAHQCRRRL

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: This antibody detects human KLHL25 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619