



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ZNF7 Polyclonal Antibody, Rabbit

Product Code: KD0101GNPAF

Gene: ZNF7 (Zinc finger protein 7, Zinc finger protein KOX4, Zinc finger protein HF.16)

Immunogen: GX5132 (GST-fusion protein, 170 amino acids)
APRTSKTDSTIRTENEQACEDMDILKSESYGTVVRISPQDFPQNPGFG
DVSDSEVWLDShLGSPGLKVTGFTFQNNCLNEETVVPKTFKTDAPQG
CKELGSSGLDCQPLESQGESAEQMSQRCEECGKGIRATSDIALHWEIN
TQKISRCQECQKKLSDCLQGKHTNNCH

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human ZNF7 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619