



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human SALL1 Polyclonal Antibody, Rabbit

Product Code: KD0111GNPAF

Gene: SALL1 (Sal-like protein 1, Zinc finger protein SALL1, Spalt-like transcription factor 1, HSal1)

Immunogen: GX5136 (GST-fusion protein, 168 amino acids)
ERTHTGEKPFACITCGRAFTTKGNLKVHMGTHMWNSTPARRGRRLSV
DGPMTFLGGNPVKFPEMFQKDLAARSGSGDPSSFWNQYAAALSNGL
AMKANEISVIQNGGIPPIPGSLGSGNSSPISGLTGNLERLQNSEPNAPLA
GLEKMASSENGTNFRFTRFVEDSKE

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human SALL1 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619