



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ILF3 Polyclonal Antibody, Rabbit

Product Code: KD0130GNPAF

Gene: ILF3 (Interleukin enhancer-binding factor 3, Nuclear factor of activated T-cells 90 kDa, NF-AT-90, Double-stranded RNA-binding protein 76, DRBP76, Translational control protein 80, TCP80, Nuclear factor associated with dsRNA, NFAR, M-phase phosphoprotein 4, MPP4)

Immunogen: GX5147 (GST-fusion protein, 170 amino acids)
DDRHVMAKHSSVYPTQEELEAVQNMVSHTERALKAVSDWIDEQEKGS
SEQAESDNMDVPPEDDSKEGAGEQKTEHMTRTLGRVMRVGLVAKGLL
LKGDLDELVLCKEKPTTALLDKVADNLAIQLAAVTEDKYEILQSVDDAA
IVIKNTKEPPLSLTIHLTSPVVREE

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human ILF3 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619