



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human TRIM69 Polyclonal Antibody, Rabbit

Product Code: KD0153GNPAF

Gene: TRIM69 (Tripartite motif-containing protein 69, RING finger protein 36, RFP-like domain-containing protein trimless)

Immunogen: GX5160 (GST-fusion protein, 166 amino acids)
KLHQFLHSKEKDILTELREEGKALNEEMELNLSQLQEQLLAKDMLVSI
QAKTEQQNSFDLKDITLLHSLEQGMKVLATRELISRKLNLGQYKGPI
QYMWREMQDTLCPGLSPLTLDPKTAHPNLVLSKSQTSVWHGDIKKIM
PDDPERFDSSVAVLGSRGFT

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human TRIM69 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619