



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ZNF192 Polyclonal Antibody, Rabbit

Product Code: KD0164GNPAF

Gene: ZNF192 (Zinc finger protein 192, LD5-1, Zinc finger protein with KRAB and SCAN domains 8)

Immunogen: GX5168 (GST-fusion protein, 100 amino acids)
KEVSARVHGHRVLWEEVVHSASAPEPPNTQLQSEATQHKSPVPQESQ
ERAMSTSQSPTRSQKGSSGDQEMTATLLTAGFQTLEKIEDMAVSLIREE
WLLD

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human ZNF192 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619