



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human ELF2 Polyclonal Antibody, Rabbit

Product Code: KD0196GNPAF

Gene: ELF2 (ETS-related transcription factor Elf-2, E74-like factor 2, New ETS-related factor)

Immunogen: GX5188 (GST-fusion protein, 152 amino acids)
VMPASTENGDKITMQPAKIITIPATQLAQCQLQTKSNLTGSGSINIVGTPL
AVRALTPVSI AHGTPVMRLSMPTQQASGQTPPRVISAVIKGPEVKSEAV
AKKQEHDVKTLQLVEEKPADGNKTVTHVVVVVSAPSAIALPVTMKTEGLV
TCE

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human ELF2 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619