



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Human COPS2 Polyclonal Antibody, Rabbit

Product Code: KD0379GNPAF

Gene: COPS2 (COP9 signalosome complex subunit 2, Signalosome subunit 2, SGN2, JAB1-containing signalosome subunit 2, Thyroid receptor-interacting protein 15, TRIP-15, Alien homolog)

Immunogen: GX5354 (GST-fusion protein, 157 amino acids)
GSPRRTTCLKYLVLANMLMKSGINPFDSQEAKPYKNDPEILAMTNLVSA
YQNNDITEFEKILKTNHSNIMDDPFIREHIEELLRNIRTQVLIKLIKPYTRIH
IPFISKELNIDVADVESLLVQCILDNTIHGRIDQVNQLLELDHQKRGGARY
TALD

Format: Affinity-purified rabbit IgG solubilized in PBS containing with 40% glycerol and 0.02% of NaN₃.

Quantity: 50 µg (250 µL)

Specificity: This antibody detects human COPS2 protein. Other species have not been tested.

Application: Immunoblotting 1:1000

Storage: -20°C

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619