



COSMO BIO CO., LTD.  
Inspiration for Life Science

## Anti-Mouse SGSM3 (FLJ00006) Polyclonal Antibody, Rabbit

**Code # :** MFL0006 [FSR03383], (Lot# XXX)

**Gene :** mouse RUN and TBC1 domain containing 3 (mSGSM3, mFLJ00006)

**Immunogen :** GX1195 (GST-fusion protein, 237 amino acids)  
FERHDDDELGFRKNDIITISQKDEHCWVGELNGLRGWFPKAFVEVLDERKEYSIAGDDSVTE  
GVTDLVRGTLCPALKALFEHGLKKPSLLGGACHPWLFIEEAAGREVERDFDSVYSRLVLCCTYR  
LDEDGKVLTPPELLYRAVQSVNVTHDAAHAQMDVKLRSLICVGLNEQVLHLWLEVLCSLPTVE  
KWYQPWSFLRSPGWVQIKCELRVLCFAFSLSQDWELPARREEEK

**Format :** Rabbit Whole Serum

**Preservative :** 0.05 % NaN<sub>3</sub>

**Quantity :** 100 µL

**Specificity :** Specific to recombinant protein GX1195. This antibody detects mSGSM3 protein. It also recognizes human SGSM3 protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at -20°C. Avoid freeze-thaw cycles.

**References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619