



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse DHX34 (KIAA0134) Polyclonal Antibody, Rabbit

Mouse KIAA0134 (mKIAA0134) protein is a homologue of human KIAA0134 (ref. 1). KIAA0134 is identical to DHX34. Rabbit anti-mouse DHX34 (mDHX34, mKIAA0134) antibody is raised against GST-fused recombinant protein (GX0622) containing following sequence.

130aa(3316-3705)

GPQTITTAPSLPGLFGNSTLSPHPTKGGYAVSDYLTYNCLTSDTDLYSDCLRSFWTCPHCGLHMPF
TPLERIAHENTCPEAPGDDPGSEEAAPPQKTSALQRPYHCQVCGQDFLFTPTEVLRHRRQHV

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0622. This antibody detects endogenous mDHX34 protein in several cell types. It also recognizes human DHX34 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>)

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MK01340910

References:

1) Okazaki, N. et al.: DNA Res., 10(4), 167 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619