



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse PDCD11 (KIAA0185) Polyclonal Antibody, Rabbit

Mouse KIAA0185 (mKIAA0185) protein is a homologue of human KIAA0185 (ref. 1). KIAA0185 is identical to PDCD11. Rabbit anti-mouse PDCD11 (mPDCD11, mKIAA0185) antibody is raised against GST-fused recombinant protein (GX0685) containing following sequence.

170aa(5188-5697)
TKSEKYKEAGELYNRMMLKRFRQEKA VVIKYGAFVLGRSQAGASHRVLQRALECLPAKEHVDVIVKFA
QLEF
QLGDVERAKAIFENTLSTYPKRTDVWSVYIDMTIKHGSQTAVRDIFERVIHLSLAPKRMKFFFKRYLDY
EKQH GTEKDVQAVKAKALEYVEAKSSALED

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0685. This antibody detects mPDCD11 protein. It also recognizes human PDCD11 protein. Details are shown in the InGaP database. (<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MK01850910

References:

1) Okazaki, N. et al.: DNA Res., 10(4), 167 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. a

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619