



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse ASTN1 (KIAA0289) Polyclonal Antibody, Rabbit

Mouse KIAA0289 (mKIAA0289) protein is a homologue of human KIAA0289 (ref. 1). KIAA0289 is identical to ASTN1. Rabbit anti-mouse ASTN1 (mASTN1, mKIAA0289) antibody is raised against GST-fused recombinant protein (GX0437) containing following sequence.

184aa(2487-3038)
LFNVQDINSALYRKVKLLNQVRLLRVQLYHMKNMFKTCRLAKELLDSFDVVPGLHTEDLHLYSLSDL
TATKKGELGPRLAELTRAGAAHVRCMLCQAKGFICEFCQNEEDVIFPFELHKCRTCEECKACYHK
TCFKSGRCPR CERLQARRELLAKQSLESYLSYDEEPEALALEATV LETT

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0437. This antibody detects mASTN1 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MK02890910

References:

1) Okazaki, N. et al.: DNA Res., 9(5), 179 (2002).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619