



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse KIAA0336 Polyclonal Antibody, Rabbit

Mouse KIAA0336 (mKIAA0336) protein is a homologue of human KIAA0336 (ref. 1, 2). KIAA0336 is identical to GRIP and coiled-coil domain-containing protein 2 (ref. 3). Rabbit anti mKIAA0336 antibody is raised against GST-fused recombinant protein (GX0491) containing following sequence.

(243 amino acids)

KVRVHNVLKQQKNKSVSQVETEGAKQEREHLEMLIDQLKIKLQDSQNSLQISVSEYQTLQAEHD
TLLERHNRMLQETVTKEAELREKLCSVQSENTMMKSEHSQTMCLTSQNEALRTSFRDQVRHL
QDEHRKTVELTQHQLSKLEAQLFQLKSEPSTRSPASSHQPSKSLRERRTTDLPLLDMMHTVAREE
GEGMETTDESSEVSSAGTHIQSLEQLLSSPDTKLERLAETSLWHNEFTKEELA

Use in any such application is responsibility of the user.

Preparation: Purified by Protein A affinity chromatography from rabbit antisera and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0491. This antibody detects endogenous levels of mKIAA0336 protein in several cell types. It also recognizes human KIAA0336 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>)

Application: Western blotting (1 : 1,000), Immunoprecipitation, Immunohistochemistry.

Storage: Keep at -20°C

Package: 100 μg (200 μL)

Product Code: MK03360310

References:

- 1) Hara, Y. *et al.*: *DNA Res.*, 10, 129 (2003).
- 2) Koga, H. *et al.*: *DNA Res.*, 11, 293 (2004).
- 3) Luke, M.R. *et al.*: *J Biol Chem.*, 278, 4216 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619