



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse PDE4DIP (KIAA0454) Polyclonal Antibody, Rabbit

Mouse KIAA0454 (mKIAA0454) protein is a homologue of human KIAA0454 (ref. 1). KIAA0454 is identical to PDE4DIP. Rabbit anti-mouse PDE4DIP (mPDE4DIP, mKIAA0454) antibody is raised against GST-fused recombinant protein (GX0271) containing following sequence.

287aa(2497-3357)

TMQSMESLLRARGLEVEQLTATCQNLQWLKEELETKFGHWQKEQESIIQQLQTS LHDRNKEVEDLSA
TLLCKLGPGQSEVAEELCQRLQRKERMLQDLLSDRNKQAVEHEMEIQGLLQSMGTREQERQAAAEK
MVQAFMERNSELQALRQYLGGKELMTSSQTFISNQPAGVTSIGPHHGEQTDQGSMQMPSRDDSTSL
TAREEASIPRSTLGDSDTVAGLEKELSSNAKEEELMAKKERESQMELSALQSMAMQEEELQVQAAD
LESLTRNVQIKEDLIK DLQMQ

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0271. This antibody detects mPDE4DIP protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MK04540910

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003)

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619