



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse-ZFP292 (KIAA0530) Polyclonal Antibody, Rabbit

Mouse KIAA0530 (mKIAA0530) protein is a homologue of human KIAA0530 (ref. 1). KIAA0530 is identical to ZFP292. Rabbit anti-mouse ZFP292 (mZFP292, mKIAA0530) antibody is raised against GST-fused recombinant protein (GX0006) containing following sequence.

136aa(5288-5695)

PMGFEASFLKFLEESAVKQKKNSDRDHSNSGSKRGSHSSRRHVDKAAVAGSSSHVCCKDSEIFV
QFANPSKLQCSENVKIVLDKTLKDRSELVLKQLQEMKPTVSLKKLEVLSNNPDRTLKEISIGKATGR
GQY

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0006. This antibody detects mZFP292 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MK05300910

References:

1) Chavrier, P. et al.: Mol. Cell. Biol. 8 (3), 1319 (1988)

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619