



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse KIAA0654 Polyclonal Antibody, Rabbit

Mouse KIAA0654 (mKIAA0654) protein is a homologue of human KIAA0654 (ref 1, 2). KIAA0654 is identical to liprin-alpha 3 (ref. 3). Rabbit anti mKIAA0654 antibody is raised against GST-fused recombinant protein (GX1885) containing C-terminal portion of mKIAA0654.

Use in any such application is responsibility of the user.

Preparation: Purified by Protein A affinity chromatography from rabbit antisera and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX1885. This antibody detects endogenous levels of mKIAA0654 protein in several cell types.

Application: Western blotting (1 : 1,000), Immunoprecipitation.

Storage: Keep at -20°C

Package: 100 µg (200 µL)

Product Code: MK06540310

References:

- 1) Hara, Y. *et al.*: *DNA Res.*, 10, 129 (2003).
- 2) Koga, H. *et al.*: *DNA Res.*, 11, 293 (2004).
- 3) Serra-Pages, C. *et al.*: *J Biol Chem.*, 273, 15611 (1998).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619