



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse FRYL (KIAA0826) Polyclonal Antibody, Rabbit

Mouse KIAA0826 (mKIAA0826) protein is a homologue of human KIAA0826 (ref. 1). KIAA0826 is identical to FRYL. Rabbit anti-mouse FRYL (mFRYL, mKIAA0826) antibody is raised against GST-fused recombinant protein (GX0380) containing following sequence.

241aa(3721-4443)

MACGLLETLKFGVLELQEHLDTYTTKREAAEQWLDNCKRTFGANEDIYRMNTNAHELEFCRRRLY
RLHFQLLLLFQAYCKLINQVNTIKNEAEVINMSEELAQLEGILKEAAASENEEIDISKAAQTTIETAI
HSLIETLKNKEFVSAVAQVKAFRTLWPNDIFGSCDDDPVQTLLHIYFHHQTLGQTGSFAVISSNLD
MSEANCKLMELNLEIRESLRTV QSYPLLAQTKPVGNMTSTGF

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0380. This antibody detects mFRYL protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MK08260910

References:

1) Okazaki, N. et al.: DNA Res., 9 (5), 179 (2002)

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619