



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse IVNS1ABP (KIAA0850) Polyclonal Antibody, Rabbit

Mouse KIAA0850 (mKIAA0850) protein is a homologue of human KIAA0850 (ref. 1). KIAA0850 is identical to IVNS1ABP. Rabbit anti-mouse IVNS1ABP (mIVNS1ABP, mKIAA0850) antibody is raised against GST-fused recombinant protein (GX0474) containing following sequence.

172aa(1853-2368)

SDPYGQKGLKNCDFVDPVTKSWTSCAPLNIRRHQSAVCELGGYLYIIGGAESWNCLNTVERYNPENNT
WTLIAPMNVARRGAGVAVLDGKLFVGGGFDGSHAISCVEMYDPTRNEWKMMGNMTSPRSNAGITTVG
NTIYAV GGFDGNEFLNTVEVYNPQSNEWSPTYKIFQF

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0474. This antibody detects mIVNS1ABP protein. It also recognizes human IVNS1ABP protein. Details are shown in the InGaP database. (<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MK08500910

References:

1) Okazaki, N. et al.: DNA Res., 10(4), 167 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619