



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse CSDE1 (KIAA0885) Polyclonal Antibody, Rabbit

Mouse KIAA0885 (mKIAA0885) protein is a homologue of human KIAA0885 (ref. 1). KIAA0885 is identical to CSDE1. Rabbit anti-mouse CSDE1 (mCSDE1, mKIAA0885) antibody is raised against GST-fused recombinant protein (GX0341) containing following sequence.

138aa(2346-2759)

QNAQTMAYNITPLRRATVECVKDQFGFINYEVGDSKKLFFHVKEVQDGVELQAGDEVEFSVILNQR
TGKCSACNVWRVCEGPKAVAAPRPDRLVNRLKNITLDDASAPRLMVLQRGPDNSMGFGAERK
IRQAGVID

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0341. This antibody detects mCSDE1 protein. It also recognizes human CSDE1 protein. Details are shown in the InGaP database. (<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MK08850910

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619