



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-Mouse ZFP644 (KIAA1221) Polyclonal Antibody, Rabbit

Mouse KIAA1221 (mKIAA1221) protein is a homologue of human KIAA1221 (ref. 1). KIAA1221 is identical to zinc finger protein 644 (ZFP644). Rabbit anti-mouse ZFP644 (mZFP644, mKIAA1221) antibody is raised against GST-fused recombinant protein (GX0214) containing following sequence:

(246 amino acids)

MMQNEEKYEKILKALNSRRIIPRFVAQKLSSGDDFLSHNVLPLDEYHNGLKTEALSVSASEEEGLH
FLSECGERKPELPSGRKNQSLTIELLKSRRLLGEERNSAVSPHKTHNQTKARKRFVQKCVLPLNEDS
PLIYQPQKMDLTMHSAIDCKQKKSRSRSGSKKMLTLPHGADDEVYILRCRFCGLVFRGPLSVQEDW
IKHLQRHIVNANLPRPGA GMVEVTSLLKKPASITETSFSLLMAEAAAS

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0214. This antibody detects mZFP644 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MK12210910

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003)

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619