



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse KIAA1250 Polyclonal Antibody, Rabbit

Mouse KIAA1250 (mKIAA1250) protein is a homologue of human KIAA1250 (ref 1, 2). KIAA1250 is identical to kinase D-interacting substance of 220 kDa (ref. 3). Rabbit anti mKIAA1250 antibody is raised against GST-fused recombinant protein (GX0506) containing following sequence.

(187 amino acids)

SDSGVRSNESSPNHSLHNEAADDSQLKANLIELEDEGHSGKRGMPHLSGLQDPVIARMSICS
EDKKSPSECSLIASSPEESWPSCQKAYNLNRTPTSTVTLNNNTAPTNRANQNFDEIEGVRETSQVI
LRPGSPNPNTAVQNENLKSMARKRSQRSSYTRLKDASELHAASSDSTGFGEERESIL

Use in any such application is responsibility of the user.

Preparation: Purified by Protein A affinity chromatography from rabbit antisera and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0506. This antibody detects endogenous levels of mKIAA1250 protein in several cell types. Details are shown in the InGaP database.

(<http://webcreate.kazusa.or.jp/create/index.jsp>)

Application: Western blotting (1 : 1,000), Immunoprecipitation.

Storage: Keep at -20°C

Package: 100 μg (200 μL)

Product Code: MK12500310

References:

- 1) Hara, Y. *et al.*: *DNA Res.*, 10, 129 (2003).
- 2) Koga, H. *et al.*: *DNA Res.*, 11, 293 (2004).
- 3) Iglesias, T. *et al.*: *J Biol Chem.*, 275, 40048 (2000).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619