



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse KBTBD2 (KIAA1489) Polyclonal Antibody, Rabbit

Mouse KIAA1489 (mKIAA1489) protein is a homologue of human KIAA1489 (ref. 1). KIAA1489 is identical to Kelch repeat and BTB (POZ) domain containing 2 (KBTBD2). Rabbit anti-mouse KBTBD2 (mKBTBD2, mKIAA1489) antibody is raised against GST-fused recombinant protein (GX0633) containing following sequence:

(233 amino acids)

GGELNRRTVERYDTEKDEWTMVSP LPCAWQWSAAVVVHDCIYVMTNLNMYCYFPRSDSWVEMAM
RQTSRSFASAAAFGDKIFYIGGLHIATNSGIRLPSGTVDGSSVTVEIYDVNKNEWKMAANIPAKRYS
CVRRAVISNLCVFMRETHLNERAKYVTYQYDLELDRLSLRQHISERVLWDLGRDFRCTVGKLYPSC
LEESPWKPTTYLFSPDGTEEFELDGEMVALPPV

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0633. This antibody detects mKBTBD2 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MK14890910

References:

1) Okazaki, N. et al.: DNA Res., 10(4), 167 (2003)

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619