



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse P140 (KIAA1684) Polyclonal Antibody, Rabbit

Mouse KIAA1684 (mKIAA1684) protein is a homologue of human KIAA1684 (ref. 1). KIAA1684 is identical to p130Cas-associated protein/SNAP-25-interacting protein (p140CAP/SNIP, P140). Rabbit anti-mouse P140 (mP140, mKIAA1684) antibody is raised against GST-fused recombinant protein (GX1199) containing following sequence:

(193 amino acids)

DLDCASKTHPGPAPTPDHKPPKAPHGQKAAPRTEPSGRRGSDELTVPRYRTEKPSKSPPPPPP
RRSFPSSHGLTTTTRTGEVVVTSKKDSVFIKKAEESELEVQKPQVKLRRRAVSEVVRPASTPPIMAS
AIKDEDDEERIIAELESGGSSVPPMKVVTPGASRLKAAQGPAGSPDKGKHGKQRTEYMRIQAQQ
Q

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX1199. This antibody detects mP140 protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MK16840910

References:

1) Okazaki, N. et al.: DNA Res., 10(4), 167 (2003)

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619