



COSMO BIO CO., LTD.  
Inspiration for Life Science

## Anti-Mouse SF3B3 (KIAA0017) Polyclonal Antibody, Rabbit

**Code # :** MKA0017PA [FIG00365], (Lot# XXX)

**Gene :** mouse splicing factor 3b, subunit 3, 130kDa (mSF3B3, mKIAA0017)

**Immunogen :** GX1030 (GST-fusion protein, 131 amino acids)  
GQEAVLAMSSRSWLSYSYQSRFHLTPLSYETLEFASGFASEQCEGIVAISTNTRLRLALEKLGA  
VFNQVAFPLQYTPRKFSVHPESNNLIIETDHNAYTEATKAQRKQQCVGSPTLVKIGMCWWVLPR  
T

**Format :** Rabbit IgG purified with Protein A affinity chromatography.

**Constitution :** PBS containing with 50% glycerol and 0.02% of NaN<sub>3</sub>

**Quantity :** 100 µg (200 µL)

**Specificity :** Specific to recombinant protein GX1030. This antibody detects mSF3B3 protein. It also recognizes human SF3B3 protein. Other species have not been tested.

**Application :** Western blotting (1 : 1,000), Other applications have not been tested.

**Storage :** Store at -20°C. Avoid freeze-thaw cycles.

**References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619