



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse CLTC (KIAA0034) Polyclonal Antibody, Rabbit

Mouse KIAA0034 (mKIAA0034) protein is a homologue of human KIAA0034 (ref. 1). KIAA0034 is identical to Clathrin heavy chain 1 (CLTC). Rabbit anti-mouse CLTC (mCLTC, mKIAA0034) antibody is raised against GST-fused recombinant protein (GX0111) containing following sequence:

(209 amino acids)

EHIGNLDRAYEFAERCNEPAVWSQLAKAQLQKGMVKEAIDSYIKADDPSSYMEVVQAANASGNW
EELVKYLQMARKKARESIVETELIFALAKTNRLAELEEFINGPNNAHIQVGDRCYDEKMYDAAKL
LYNNVSNFGRLASTLVHLGEYQAAVDGARKANSTRTWKEVCFACVDGKEFRLAQMCGHLHVVHA
DELEELINYYQDRGF

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0111. This antibody detects mCLTC protein. Details are shown in the InGaP database.
(<http://webcreate.kazusa.or.jp/create/index.jsp>) Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MKA0034AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619