



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse SETDB1 (KIAA0067) Polyclonal Antibody, Rabbit

Mouse KIAA0067 (mKIAA0067) protein is a homologue of human KIAA0067 (ref. 1). KIAA0067 is identical to SET domain, bifurcated 1 (SETDB1). Rabbit anti-mouse SETDB1 (mSETDB1, mKIAA0067) antibody is raised against GST-fused recombinant protein (GX0398) containing following sequence:

(156 amino acids)

ISSGSDGDDDFEDKKNLSGPTKRQVAVKSTRGFALKSTHGIAIKSTNMAVDKGESAPVRKNTRQFYD
GEESCYIIDAKLEGNLGRYLNHSCSPNLFVQNVFVDTHDLRFPWVAFFASKRIRAGTELTWDYNYEV
GSVEGKELL CCCGAIECRGRL

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0398. This antibody detects mSETDB1 protein. It also recognizes human SETDB1 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MKA0067AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619