



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse SERTAD2 (KIAA0127) Polyclonal Antibody, Rabbit

- Code # :** MKA0127 [FSR00419], (Lot# XXX)
- Gene :** mouse SERTA domain containing 2 (mSERTAD2, mKIAA0127)
- Immunogen :** GX2275 (GST-fusion protein, 120 amino acids)
AAHTAAPEGPKGTSSSESVQKPEGPEEGRTDDSRFMDSLPGNFEITTSTGFLDTLDDILFADI
DTSMYDFDPCTSASGTASKMAPVSADDLLKTLAPYSNQPVAPSQPFKMDLTLDH
- Format :** Rabbit Whole Serum
- Preservative :** 0.05 % NaN₃
- Quantity :** 100 µL
- Specificity :** Specific to recombinant protein GX2275. This antibody detects mSERTAD2 protein. It also recognizes human SERTAD2 protein. Other species have not been tested.
- Application :** Western blotting (1 : 1,000), Other applications have not been tested.
- Storage :** Store at -20°C. Avoid freeze-thaw cycles.
- References :** 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by  ProteinExpress Co., Ltd. and  Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619