



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse TTLL12 (KIAA0153) Polyclonal Antibody, Rabbit

Code # : MKA0153 [FSR00432], (Lot# XXX)

Gene : mouse tubulin tyrosine ligase-like family, member 12 (mTTLL12, mKIAA0153)

Immunogen : GX2045 (GST-fusion protein, 151 amino acids)
RFSNRPFALDDLDYKHTVMNYDPDVVLKQVHYNEFIPQFEKQYPEFPWSDVQAEIFKAFTE
LFQVACAKPPPMGLCDYPSSRAMYAIDLMLNWDNHPDGKRVMQPQILEVNFNPDCERACRYH
PSFFNDVFSTLFLDETNDCHVTRII

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX2045. This antibody detects mTTLL12 protein. It also recognizes human TTLL12 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619