



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse SCR1 (KIAA0193) Polyclonal Antibody, Rabbit

Code # : MKA0193 [FSR00450], (Lot# XXX)

Gene : mouse secernin 1 (mSCR1, mKIAA0193)

Immunogen : GX0876 (GST-fusion protein, 224 amino acids)
LSLATKLDEEHPELRTYAQSQGWWTGDDEFNFAQVFPADRLDCCAGQDSLEKQEEESITVQ
TMINILRDKASGVCIDSESFLTASIVSVLPQNRSSPCIHFTGTPDPSRSIFKPFIFVDDVKLVPK
AQSPCFGDDDDPAKKEPRFQEKPDRRHELYKAHEWARAVIESDEEQGRTLRLKTMLELEKQGLE
AMDEILSSPEPPDPAEVGDLFYDCVDTEMKFFK

Format : Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity : 100 µL

Specificity : Specific to recombinant protein GX0876. This antibody detects mSCR1 protein. It also recognizes human SCR1 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619