



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse SCRN1 (KIAA0193) Polyclonal Antibody, Rabbit

Code # : MKA0193PA [FIG00450], (Lot# XXX)

Gene : mouse secernin 1 (mSCRN1, mKIAA0193)

Immunogen : GX0876 (GST-fusion protein, 224 amino acids)
LSLATKLDEEHPELRTYAQSQGWWTGDDEFNFAQVFPADRLDCCAGQDSLEKQEEESITVQ
TMINILRDKASGVCIDSESFLTASIVSVLPQNRSSPCIHFTGTDPSPRSIFKPFIFVDDVKLVPK
AQSPCFGDDDDPAKKEPRFQEKPDRRHELYKAHEWARAVIESDEEQGRTLRLKTMLELEKQGLE
AMDEILSSPEPPDPAEVGDLFYDCVDTEMKFFK

Format : Rabbit IgG purified with Protein A affinity chromatography.

Constitution : PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity : 100 µg (200 µL)

Specificity : Specific to recombinant protein GX0876. This antibody detects mSCRN1 protein. It also recognizes human SCRN1 protein. Other species have not been tested.

Application : Western blotting (1 : 1,000), Other applications have not been tested.

Storage : Store at -20°C. Avoid freeze-thaw cycles.

References : 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by  ProteinExpress Co., Ltd. and  Kazusa DNA Research Institute.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618 FAX : +81-3-5632-9619