



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse PHF16 (KIAA0215) Polyclonal Antibody, Rabbit

Mouse KIAA0251 (mKIAA0251) protein is a homologue of human KIAA0251 (ref. 1). KIAA0251 is identical to PHD finger protein 16 (PHF16). Rabbit anti-mouse PHF16 (mPHF16, mKIAA0251) antibody is raised against GST-fused recombinant protein (GX1028) containing following sequence:

(162 amino acids)

RVSSSNGLEGNWSGNITQKVNSSSEVCYDQESMLSSHLPSPGNIRKSSMEHFSSRSFKEATNTWV
KPTEDLQYCVKPTKNVSSKEQLWGRQLLRPTGRASYQETDGYCPDLEPSDSEAEGEGSKETP
RVKRESSDRENPS HDSARECHGKTKTHPHSHSSMQR

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX1028. This antibody detects mPHF16 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000).

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA0215AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use



ProteinExpress Co., Ltd. and

Manufactured by



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619