



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-Mouse RIMS3 (KIAA0237) Polyclonal Antibody, Rabbit

Product Code: PRX-MKA0237

Gene: mouse regulating synaptic membrane exocytosis 3 (mRIMS3, mKIAA0237)

Immunogen: GX1817 (GST-fusion protein, 115 amino acids)

YIKAYLLENGACVAKKTKVAKKTCDPYQQALLFDEGPQGKVLQVIVWGDYGRMDHKCFMG

MAQIMLDELDSLAAVTGWYKLFPTSSVADSTLGSLTRRLSQSSLESATSPSCS

ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA0237](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: Specific to recombinant protein GX1817. This antibody detects mRIMS3 protein. It also recognizes human RIMS3 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619