



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti- Mouse CCDC132 (KIAA1861) Polyclonal Antibody, Rabbit

Product Code: PRX-MKA0313PA

Gene: mouse Rap guanine nucleotide exchange factor (GEF) 2
(mRAPGEF2, mKIAA0313)

Immunogen: GX0190 (GST-fusion protein, 107 amino acids)
NNLEREKMGGHLRLLNIACAIAKAKRRRLMTLTKPSREAPLPFILLGGSEKGFIFVDSVDSCKAT
EAGLKRGDQILEVNGQNFENIQLSKAMEILRNNTLSITVKT

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit IgG purified with Protein A affinity chromatography.

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 100 µg (200 µL)

Specificity: Specific to recombinant protein GX0190. This antibody detects mRAPGEF2 protein. It also recognizes human RAPGEF2 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619