



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-Mouse CLEC16A (KIAA0350) Polyclonal Antibody, Rabbit

Product Code: PRX-MKA0350PA

Gene: mouse C-type lectin domain family 16, member A (mCLEC16A, mKIAA0350)

Immunogen: GX1281 (GST-fusion protein, 302 amino acids)

QNTTDEEKSAATNSENAQWSRPFLDMVYHALDSPDDDYHALFVLCCLYAMSHNKGMDPEKLK
RIQLPVPSEAEKTTYNHLLAERLIRIMNNAQPDGRIRLATLELSCLLLKQQVLTSSGCVIKDVHL
ACLEGAREESVHLVRHFYKGEEIFLDMFEDEYRSMTIKPMNVEYLMMDASILLPPTGTPLTGIDF
VKRLPCGDVEKTRRAIRVFFMLRSLQLRGEPETQLPLTREEDLIKTDDVLDLNNSDLIACTVIT
KDGGMVQRFLAVDIYQMSLVEPDVSRLLGWGVVKFAGLLQDMQV

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit IgG purified with Protein A affinity chromatography.

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 100 µg (200 µL)

Specificity: Specific to recombinant protein GX1281. This antibody detects mCLEC16A protein. It also recognizes human CLEC16A protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619