



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti- Mouse MADD (KIAA0358) Polyclonal Antibody, Rabbit

Product Code: PRX-MKA0358

Gene: mouse MAP kinase-activating death domain protein, Differentially expressed in normal and neoplastic cells, Insulinoma glucagonoma clone 20, Rab3 GDP/GTP exchange factor (mMADD, mKIAA0358)

Immunogen: GX0189 (GST-fusion protein, 252 amino acids)

SPAVMIRSSSQDSEVSTVSNSSGETLGADSDLSSNAGDGPGGEGSAHLASSRATLSDSEIETN
SATSAIFGKAHSLKPKEKPAGSPIRSEDVSQRVVLYEGLLGRDKGSMWDQLEDAAMETFSLS
KERSTLWDQMWFEDAFDAVMLEREGMGMDQGPQEMIDRYLSLGEHDRKRLEDDERLLA
TLLHNLISYMLLMKVNKNDIRKKVRRMLMGKSHIGLVYSQQVNEVLDQLNSLNGRDLISRSSGSR
H

ROUGE Database (Mouse KIAA/FLJ cDNA Database)

[mKIAA0358](#)

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit Whole Serum

Preservative : 0.05 % NaN₃

Quantity: 100 µL

Specificity: Specific to recombinant protein GX0189. This antibody detects mMADD protein. It also recognizes human MADD protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619