



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse KIF3B (KIAA0359) Polyclonal Antibody, Rabbit

Mouse KIAA0359 (mKIAA0359) protein is a homologue of human KIAA0359 (ref. 1). KIAA0359 is identical to kinesin family member 3B (KIF3B). Rabbit anti-mouse KIF3B (mKIF3B, mKIAA0359) antibody is raised against GST-fused recombinant protein (GX2286) containing following sequence:

(158 amino acids)

HSLVAEEKMRLLKEKEKKMEDLRREKDAAEMLGAKIKAMESKLLVGGKNIVDHTNEQQKILEQKR
QEIAEQKRREREIQQMESRDEETLELKETYTSLQQEVDIKTKLKLFSKLQAVKAEIHDLQEEHI
KERQELEQTQNEL TRELKLKHLIEN

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX2286. This antibody detects mKIF3B protein. It also recognizes human KIF3B protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Immunoprecipitation.

Storage: Keep at -20°C

Package: 50 µg (250 µL)

Product Code: MKA0359AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619