



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse MYST4 (KIAA0383) Polyclonal Antibody, Rabbit

Mouse KIAA0383 (mKIAA0383) protein is a homologue of human KIAA0383 (ref. 1). KIAA0383 is identical to MYST histone acetyltransferase (monocytic leukemia) 4 (MYST4). Rabbit anti-mouse MYST4 (mMYST4, mKIAA0383) antibody is raised against GST-fused recombinant protein (GX1056) containing following sequence:

(159 amino acids)

QLELSVQDGSVLKVTNKGLASYKDPDNPGRFSSVKPGTFPKPTKGSKGPPCNDLRNVDWNKLLK
RAIEGLEEPNGSSSLKNIEKYLRSQSDLTGTTNHPAFQQRLRLGAKRAVNNGRLLKEGPQYRVNSG
SSDGKGAPQYPS AFPSSLPPVSLLPHEKDQ

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX1056. This antibody detects mMYST4 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA0383AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619