



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-Mouse MYO6 (KIAA0389) Polyclonal Antibody, Rabbit

Product Code: PRX-MKA0389PA

Gene: mouse myosin VI (mMYO6, mKIAA0389)

Immunogen: GX1161 (GST-fusion protein, 163 amino acids)

NTETEQRAPKSVTDYDFAPFLNNSPQQNPAAQLPARQQEIDMKRQQRFFRIPFIRPADQYKDP
QNKKKGWWYAHFDGPWIARQMELHPDKPPILLVAGKDDMEMCELNLEETGLTRKRGAEILPR
QFEEIWERCGGIQYLQSAIESRQARPTYATAMLQNLLK

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit IgG purified with Protein A affinity chromatography.

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 100 µg (200 µL)

Specificity: Specific to recombinant protein GX1161. This antibody detects mMYO6 protein. It also recognizes human MYO6 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619