



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-Mouse DDEF2 (KIAA0400) Polyclonal Antibody, Rabbit

Mouse KIAA0400 (mKIAA0400) protein is a homologue of human KIAA0400 (ref. 1). KIAA0400 is identical to development and differentiation enhancing factor 2 (DDEF2). Rabbit anti-mouse DDEF2 (mDDEF2, mKIAA0400) antibody is raised against GST-fused recombinant protein (GX1332) containing following sequence:

(251 amino acids)

YSRMQSLTLDVLGTSELLAKNIGNAGFNEIMECCLPSEDPVKPNPGSDMIARKDYITAKYMERRY
ARKKHADTAAKLHSLCEAVKTRDIFGLLQAYADGVDLTEKIPLANGHEPDETALHLAVRSVDRTSLHI
VDFLVQNSGNLDKQGTGKGSTALHYCCCLTDNAECLKLLLRGKASIEIANESGETPLDIARKLKHEHCE
ELLTQALSGRFNSHVHV EYEWRLLEDLDESDDVDKLPSPNRREDRP

This product is for laboratory use only; use in any such application is responsibility of the user.

Preparation: Affinity-purified from rabbit antiserum and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX1332. This antibody detects mDDEF2 protein. It also recognizes human DDEF2 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000)

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MKA0400AF

References:

1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Res.Inst.

Kazusa DNA Research Institute.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619