



COSMO BIO Co., LTD.
Inspiration for Life Science

Anti-Mouse NUPL1 (KIAA0410) Polyclonal Antibody, Rabbit

Product Code: PRX-MKA0410PA

Gene: mouse nucleoporin like 1 (mNUPL1, mKIAA0410)

Immunogen: GX0924 (GST-fusion protein, 269 amino acids)

LKVQEDIKALKQLLSLAASGLQRNTLNIDKLETAQELKNAEIALRTQKTPPGLQHENTAPADYF
RILVQQFEVQLQQYRQQIEELENHLATQASNSHITPDLSMAMQKIYQTFVALAAQLQSIHENVK
VLKEQYLG YRKMFLGDAVDVFEARRTEAKKWQNAPRVTTGPTPFSTMPNAAVAMAATLTQQ
QQPATGPQPSLGVSGFTPFSGIGTGLQSSGLSSNLGGFGTSSGFGCGTTGASTFGFGTTD
KPSGSLSAGLCASA

The mammalian expression vector of the full length human protein corresponding to the antibody is available!

[Human cDNA/ORF Clones in Kazusa Collection \(KDRI\)](#)

Format: Rabbit IgG purified with Protein A affinity chromatography.

Constitution: PBS containing with 50% glycerol and 0.02% of NaN₃

Quantity: 100 µg (200 µL)

Specificity: Specific to recombinant protein GX0924. This antibody detects mNUPL1 protein. It also recognizes human NUPL1 protein. Other species have not been tested.

Application: Western blotting (1 : 1,000), Other applications have not been tested.

Storage: Store at -20°C. Avoid freeze-thaw cycles

References 1) Okazaki, N. et al.: DNA Res., 10(1), 35 (2003).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619